









Active Electrode

Ball Electrode















Ø2.38mm

Description	Package	Ref. No.
*Ball Ø4.0mm * 50mm non-stick electrode with shrinking and molding, sterile	 =100	SW130503121S SW130503111S (bulk)
		Scale 1:1
*Ball Ø5.0mm * 70mm non-stick electrode with shrinking, bulk	 =100	SW130703118S
*Ball Ø5.0mm * 70mm non-stick electrode with shrinking, sterile	 = 50	SW130703119S
		Scale 1:1
*Ball Ø5.0mm * 70mm non-stick electrode with shrinking & molding, sterile	 = 50	SW130703139S
*Ball Ø5.0mm * 70mm non-stick electrode with shrinking & molding, bulk	 =100	SW130703138S
		Scale 1:1
Ball Ø5.0mm * 70mm electrode with shrinking, bulk	 =100	SW130703217S
Ball Ø5.0mm * 70mm electrode with shrinking, sterile	 = 50	SW130703218S
		Scale 1:1
Ball Ø5.0mm * 70mm electrode with shrinking & molding, sterile	 = 50	SW130703242S
		Scale 1:1
*Ball Ø3.0 * 110mm non-stick electrode with shrinking & molding, bulk	 = 50	SW131103111S
		Scale 1:1

Active Electrode

Ball Electrode

Ø2.38mm

Description	Package	Ref. No.
*Ball Ø5.0 * 110mm non-stick electrode with shrinking & molding, bulk	 =100	SW131103112S
*Ball Ø5.0 * 110mm non-stick electrode with shrinking & molding, sterile	 = 50	SW131103131S
		
Scale 1:1		
Ball Ø4.0 * 130mm electrode with shrinking, sterile	 = 50	SW131303211S
		
Scale 1:1		
Ball Ø4.0 * 130mm electrode with shrinking & molding, sterile	 = 50	SW131303212S
		
Scale 1:1		
*Ball Ø5.0 * 150mm non-stick electrode with shrinking & molding, bulk	 =100	SW131503137S
*Ball Ø5.0 * 150mm non-stick electrode with shrinking & molding, sterile	 = 50	SW131503138S
		
Scale 1:1		
Ball Ø5.0 * 150mm electrode with shrinking, bulk	 =100	SW131503219S
Ball Ø5.0 * 150mm electrode with shrinking, sterile	 = 50	SW131503220S
		
Scale 1:1		
Ball Ø5.0 * 150mm electrode with shrinking & molding, bulk	 =100	SW131503237S
Ball Ø5.0 * 150mm electrode with shrinking & molding, sterile	 = 50	SW131503238S
		
Scale 1:1		