












電燒配件

雙極鑷子電燒線

母端 (器械)	公端 (電燒機)		
		ø 4mm x 2 prong pin	28.8mm Valleylab/Bowa/Aesculap Lamideay/Olympus/Soering
 <p>European type 3.8 6.7 European Type 3.8mm Flat Plug Length: 3M</p>		SW4643002A-T (灰)	SW4643001E (灰)
 <p>European type 4.2 6.7 European Type 4.2mm Flat Plug Length: 3M</p>		SW4643004A (藍)(滅菌) SW4643004B (藍)(大宗包)	SW4643004G (藍)(滅菌) SW4643004D (藍)(大宗包)
 <p>US type 5 2 US Type ø 2.0mm Angled Plug Length: 3M</p>		SW4613005A-T (藍)(滅菌) SW4613005E (藍)(大宗包) SW4613003E (灰)	SW4613005B (藍)(滅菌) SW4613005D (藍)(大宗包) SW4613003F (灰)
 <p>US Type ø 2.0mm Protect Plug Length: 3M</p>		SW4613003I (灰)	SW4613003H (灰)
 <p>Aesculap Type Protect Plug Length: 3M</p>		SW4643007B (藍)(大宗包) SW4643007C (灰)	SW4643007A (灰)

※ (灰)：重覆式 / (藍)：拋棄式

雙極鑷子電燒線

<p>公端 (電燒機)</p>	
<p>母端 (器械)</p>	<p>ø 4mm x 2 prong pin</p>
 <p>US Type Ø2mm Non Booted plug Length: 3M</p>	<p>SW4613003I-01(灰)</p>
 <p>Aesculap Type Non Booted Plug Length: 3M</p>	<p>SW4643007D(灰)</p>
 <p>2.0mm x 2 Flat Plug</p>	<p>SW4613004B(灰)</p>

Description

	Feature	Color	Type	Description		Length	Material	Sterile
				Male Plug	Female Plug			
				(Generator side)	(Instrument side)			
SW4613005A-T	Disposable	Blue	US	4mm x 2 prong pin	2.0mm flat plug	3 M	PVC	Y
SW4613005E	Disposable	Blue	US	4mm x 2 prong pin	2.0mm flat plug	3 M	PVC	N
SW4613005B	Disposable	Blue	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	PVC	Y
SW4613005D	Disposable	Blue	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	PVC	N
SW4613003F	Reusable	Grey	US	28.8mm 2 prong pin	2.0mm flat plug	3 M	PVC	N
SW4643004A	Disposable	Blue	European	4mm x 2 prong pin	4.2mm flat plug	3 M	PVC	Y
SW4643004B	Disposable	Blue	European	4mm x 2 prong pin	4.2mm flat plug	3 M	PVC	N
SW4643004D	Disposable	Blue	European	28.8mm 2 prong pin	4.2mm flat plug	3 M	PVC	N
SW4643004G	Disposable	Blue	European	28.8mm 2 prong pin	4.2mm flat plug	3 M	PVC	Y
SW4613003E	Reusable	Grey	US	4mm x 2 prong pin	2.0mm flat plug	3 M	TPR plug + Silicon cord	Y
SW4613003H	Reusable	Grey	US	28.8mm 2 prong pin	2.0mm protect plug	3 M	TPR plug + Silicon cord	N
SW4613003I	Reusable	Grey	US	4mm x 2 prong pin	2.0mm protect plug	3 M	TPR plug + Silicon cord	N
SW4643001E	Reusable	Grey	European	28.8mm 2 prong pin	3.8mm flat plug	3 M	TPR plug + Silicon cord	N
SW4643002A-T	Reusable	Grey	European	4mm x 2 prong pin	3.8mm flat plug	3 M	TPR plug + Silicon cord	N
SW4643007A	Reusable	Grey	Aesculap	28.8mm 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	N
SW4643007B	Reusable	Blue	Aesculap	4mm x 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	N
SW4643007C	Reusable	Grey	Aesculap	4mm x 2 prong pin	Flat plug	3 M	TPR plug + Silicon cord	N
SW4613003I-01	Reusable	Grey	US	4mm x 2 prong pin	2.0 mm non booted plug	3 M	TPR plug + Silicon cord	N

SW4643007D	Reusable	Grey	Aesculap	4mm x 2 prong pin	Flat non booted plug	3M	TPR plug + Silicon cord	N
SW4613004B	Reusable	Grey		4mm x 2 prong pin	2.0mm x 2 flat plug	3M	TPR plug + Silicon cord	N



升望科技股份有限公司 Shining World Health Care Co., Ltd.

電話：02-22900966 傳真：02-22903966 www.shinmed.com shinmed17@shinmed.com